

Generating Nanoparticles from Tellurium, Bismuth, Gallium, Gold, Iron and Various Salts Using the Silver Particle Generator (SPG)


V. Berger¹, K.M. Risby², P. S. Bauer¹, M. Dollner¹, H.J. Schulz¹, J. Swanson^{1,3}, A. Boies^{1,4}

¹ Catalytic Instruments, ² University of Cambridge, ³ Minnesota State University, ⁴ Stanford University

APT, Corfu, 11 May 2026

Aerosol Generation


- Start of many lab studies in aerosol science:




- Evaporation-condensation method
 - Scheibel and Porstendörfer (1983). (Salt (NaCl) and Silver)
 - versatile and reliable technique
 - Different materials possible: Metals, Salts, Plastics (Wlasits *et al.*, 2022)
- Each lab has a different setup → difficult to compare
- Tube furnaces cumbersome to operate, slow temperature changes

Silver Particle Generator (SPG)

- Standardize and simplify the evaporation-condensation generation method
- Precise, fast and easy control of temperature and flow
- Stable and high repeatability
- Different materials are possible
- Simple exchange of tube set




SPG Tube Set



- Ceramic (Alumina) Tubes
- Cleaned by baking
- 1 Tube set for each material


CPC Counting Efficiency for Different Materials

- The material of seed particles influence condensation behavior
- Example:
Condensation Particle Counter (CPC)
Counting efficiency
 - Krasa, H., Risby, K.M., Kupper, M., Bergmann, A., and Boies, A.M. (2025). *Aerosol Science and Technology* 59 (9):1151



Krasa et al. (2025). *Aerosol Science and Technology* 59 (9):1151



Low Temperature Materials (up to 300 °C)



- Tetracontane ($\text{C}_{40}\text{H}_{82}$)
 - PMP Volatile Particle Remover Testing
 - Testing our own Catalytic Strippers (CS)
- Plastics
(similar as Wlasits *et al.*, 2022)

Medium Temperature Salts (300 – 900 °C)


Sodium iodide (NaI)




Krasa *et al.* (2025). *Aerosol Science and Technology* 59 (9):1151

Potassium Chloride (KCl)

KCl Geometric Mean Diameter (GMD)




KCl Vapor Pressure




Potassium Chloride (KCl)



KCl Geometric Mean Diameter (GMD)





Geometric Mean Diameter (GMD)





Potassium Chloride (KCl)



Potassium Chloride (KCl)




Salts at Medium Temperatures




Krasa et al. (2025). *Aerosol Science and Technology* 59 (9):1151


Medium Temperature Materials





Medium Temperature Materials




Medium and High Temperature Materials



Silver Particles




SINTERING
S8000 @ 400 °C



≈ 100 nm


- solid, rather unreactive, normally insoluble
- can be formed to a thermally stable and compact sphere (sintering)

Medium and High Temperature Materials



Conclusion

- SPG generates aerosol of various materials by evaporation-condensation
- Precise, fast and easy control
- Stable and high repeatability
- Simple exchange of tube sets
- Contact us with your material needs!



°Catalytic Instruments®

hot technologies • clean solutions




Contact us

Elevate Your Aerosol
Measurement Precision with us!
sales@catalytic-instruments.com


Thank you for your attention!



Silver on Gold seed particles




Experimental Setup



- Gold seed particles routed into silver SPG for heterogeneous nucleation
- Tested homogeneous nucleation with HEPA filter

Tandem SPG experiments



Tandem SPG experiments

